
**bio silver**

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
```
OK | Started Remount Root and Kernel File Systems.
    Starting Load/Save Random Seed...
OK | Started Load/Save Random Seed.
    Starting Load/Save Random Seed...
OK | Started Entropy daemon using the HWRNGE algorithm.
    Starting Journal Service...
OK | Started Monitoring of LVM2 mirrors, snapshots etc. using dmeventd or progress polling.
OK | Reached target Local File Systems (Pre).
OK | Started Journal Service.
    Starting Flush Journal to Persistent Storage...
OK | Started Flush Journal to Persistent Storage...
OK | Started udev Coldplug all Devices.
OK | Started Dispatch Password Requests to Console Directory Watch.
OK | Found device SanDisk_SSD_PLUS_120_GB 3.
OK | Activating swap Swap Partition...
OK | Activated swap Swap Partition.
OK | Reached target Swap.
OK | Found device SanDisk_SSD_PLUS_120_GB EFI\x20System\x20Partition.
    Starting File System Check on /dev/disk/by-uuid/29DA-5B1A...
OK | Started File System Check Daemon to report status.
OK | Started File System Check on /dev/disk/by-uuid/29DA-5B1A.
    Mounting /boot/efi...
OK | Mounted /boot/efi.
OK | Reached target Local File Systems.
    Starting Set console font and keypad...
    Starting Tell Plymouth To Write Out Runtime Data...
    Starting Clean up any mess left by odms-up...
    Starting Create Volatile Files and Directories...
    Starting LSB: AppArmor initialization...
    Starting Run snappy Firstboot setup...
OK | Listening on Load/Save RF Kill Switch Status /dev/rfkill Watch.
OK | Started Tell Plymouth To Write Out Runtime Data.
OK | Started Create Volatile Files and Directories.
OK | Started Clean up any mess left by odms-up
```

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The login screen for smilebaby features a large 'smilebaby' logo at the top. Below it, there is a 'ログイン' (Login) section with two input fields: 'ユーザーID' (User ID) containing 'smilebaby' and 'パスワード' (Password) with masked characters. A checkbox labeled 'ログイン情報を記録する' (Remember login information) is checked. A 'ログイン' (Login) button is positioned at the bottom right of the form.

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The dashboard shows the '見守り一覧' (Monitoring List) page. At the top, there are navigation tabs for '見守り一覧', 'ベッド情報設定', and 'アラート履歴'. Below the tabs, there are search filters for '建物名', 'フロア名', and '部屋名', along with a '利用者名' field and '検索' (Search) and 'クリア' (Clear) buttons. A 'ベッド数' (Bed count) section shows '1' and a '表示フィルター' (Display filter) section with buttons for '開始', '停止', '着床', and '離床'. On the right, there are '全開始' (All Start) and '全停止' (All Stop) buttons. The main table has columns for '見守り', '着床', '部屋', '名前', '心拍', '呼吸', '体動', '離床時間', 'アラート', and 'グラフ'. The first row shows a checked '見守り' status, a bed icon, '利用者1', heart rate '0', breathing '0', and a '離床' alert at '02/25 14:04:08'. There are also icons for '離床', '呼吸', '心拍', and '体動' in the alert column.

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